

**Baxter**

Oxiris

CYTOKINE - ENDOTOXIN AND  
CRRT HEMOFILTER



## Oxiris Set

### 3-IN-1 CRRT-SEPSIS MANAGEMENT

Blood purification beyond CRRT  
by targeting cytokine and endotoxin removal Powered By  
**PrisMax** and **Prismaflex**


# OXIRIS SET

The Oxiris set is indicated for use only with the **Prismaflex** control unit or with **PrisMax** control unit (in countries where **PrisMax** is cleared or registered). It is intended for patients in need of blood purification, including continuous renal replacement therapy, and in conditions where excessive endotoxin and inflammatory mediator levels exist.

OXIRIS SET GENERAL DATA					
Weight		890 g			
Overall Dimensions		27 x 22 x 9 cm			
Blood volume in set ± 10 %		193 ml			
Minimal patient weight		30 kg			
Materials					
Oxiris hollow fiber: Acrylonitrile and sodium methallyl sulfonate copolymer + Polyethylenimine (surface treatment agent) + heparin grafted [4500+/- 1500 IU/m²]					
Filter housing and headers: Polycarbonate					
Filter potting compound: Polyurethane					
Tubing material: Plasticized polyvinyl chloride (PVC)					
Cartridge: PETG					
Sterilization mode: EtO (ethylene oxide)					
Filter operating specifications					
Maximum TMP* (mmHg/kPa)		450/60			
Maximum blood pressure (mmHg/kPa)		500/66.6			
Range of blood flow rate		100-450 ml/min			
Filter data					
Nominal physical characteristics:					
Effective surface area		1.5 m²			
Fiber internal diameter (wet)		240 µm			
Fiber wall thickness		50 µm			
IN VITRO PERFORMANCES					
CVVHD clearances					
Clearances versus inlet dialysate flow rate (Continuous veno-venous hemodialysis) (Saline, T = 37°C).					
		Oxiris set Q <sub>B</sub> ** = 200 ml/min Q <sub>UF</sub> **** = 0 ml/min			
QD	l/h ml/min	1 17	2.5 42	4 67	8 133
Urea (±10%)		17	42	66	117
Vitamin B <sub>12</sub> (±20%)		17	38	51	68
Inulin (±20%)		16	33	40	49
*Transmembrane pressure. **Access blood flow rate. ***Protein concentration. ****Ultrafiltration flow rate <sup>(1)</sup> . <sup>(1)</sup> The ultrafiltration flow rate is the "patient fluid removal flow rate + replacement flow rate + pre-blood-pump flow rate".					

## CVVH performances<sup>1</sup>

"In vitro" ultrafiltration with blood (values ±15%)  
(Continuous veno-venous hemofiltration)  
(Bovine blood at 37°C, Hematocrit 32%, Cp<sup>\*\*\*</sup> 60 g/l).  
Ultrafiltration is controlled by the **Prismaflex** system and is independent of the ultrafiltration coefficient (KUF)



## Sieving coefficient

(Bovine plasma, Cp 60 g/l, T = 37°C)  
Q<sub>B</sub> = 100 ml/min, Q<sub>UF</sub> = 20 ml/min

Urea	1
Vitamin B <sub>12</sub>	1
Inulin	0.96
(Human plasma, Cp 60 g/l, T=37°C)	
Myoglobin	0.70
Albumin	<0.0045

## Cytokine adsorption

Cytokine adsorption removal rate (%) <sup>(2)</sup>  
(human plasma, Cp 60 g/l, 37°C)  
Q<sub>B</sub> = 150 ml/min, Q<sub>UF</sub> = 0 ml/min

IL-10 (± 10%)	96
IL-6 (± 10%)	84
HMGB-1 (± 10%)	94
TNF-α (± 30%)	82

<sup>(2)</sup>Removal Rate expressed at t=120 min with a theoretical initial IL-10, IL-6, HMGB-1 and TNF-α respective concentration of 500 pg/ml, 1500 pg/ml, 30 ng/ml and 250 pg/ml.

## Endotoxin adsorption

Lipopolysaccharide adsorption removal rate (%)<sup>(3)</sup>  
(human plasma, Cp 60 g/l, 37°C)  
Q<sub>B</sub> = 150 ml/min, Q<sub>UF</sub> = 0 ml/min

LPS (± 20%)	75
-------------	----

<sup>(3)</sup>Removal Rate expressed at t=120 min with an initial LPS concentration after stabilization of 50±10 EU/ml  
Cp: Protein concentration  
RR: removal rate  
IL-10: Interleukin-10  
IL-6: Interleukin-6  
HMGB-1: High-mobility group box 1  
TNF-α: Tumor necrosis factor – α  
LPS: Lipopolysaccharide

## ORDERING INFORMATION

	Code N°	N° units/box
Oxiris S set	955503	4

1. Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions. Adsorption removal rate obtained in vitro are likely to differ from in vivo results. Adsorption characteristics change with the duration of observation.

For safe and proper use of the devices mentioned herein, please refer to the Instructions for Use.

Baxter, Oxiris, Prismaflex and Prismax are trademarks of Baxter International Inc. or its subsidiaries.  
NORD/MG146/20-0003 – August 2020

[www.baxter.se](http://www.baxter.se)  
Baxter Medical AB  
P 08 632 64 00

[www.baxter.no](http://www.baxter.no)  
Baxter AS  
P 22 58 48 00

[www.baxter.fi](http://www.baxter.fi)  
Baxter Oy  
P (09) 8621 111

[www.baxter.dk](http://www.baxter.dk)  
Baxter A/S  
P 4816 6400